



PRODUCT BULLETIN

Product(s):	COMHUB
Bulletin Number	040
Release Date:	March 17, 1995
Pages:	2

Subject: COMHUB Software Release 1.05

This COMHUB Software Release (1.05) contains two main features requested by our customers:

- (1) Datascope Multinex Monitor support
- (2) Auto configure of CO2 units

Datascope Multinex Support

The Multinex feature allows the following Datascope models to interface with all Narkomed's via connection with the COMHUB and remote display:

	SPO2	CO2 Analysis	Agent Analysis	Agent Selection
MODEL 4000	X	X	X	Manual
MODEL 4000ID	X	X	X	Auto - ID
MODEL 4100		X	X	Manual
MODEL 4100ID		X	X	Auto - ID
MODEL 4200	X	X		
MODEL 4300		X		

With these monitors the following information will be displayed on the Remote Display Screens and stored in the Data Management record: **SPO2 Sat, SPO2 Pulse, CO2, O2, Agent Data.**

Auto Configure of CO2

Currently the units of measure for CO2 are set at mmHg and must be manually matched to the Datascope monitor units of measure if different. The Auto configure of CO2 units allows the user to change units at the third party monitor and the CO2 units will be automatically changed on the NM2C and NM2B (NM4 will follow with a future software upgrade). For the NM2B, the ORDM software requires an update to support this automatic function.

MULTINEX CABLE PART NUMBER AND PRICE

The part number for the Datascope Multinex cable is #4113142. The List Price is \$85.00.

SERVICE SOFTWARE TOOLKIT (SST)

SST Software version 1.04 will be required to support the Multinex configuration.

DATASCOPE UPDATE

There is an update required for the Datascope monitor in order to support Desflurane and Sevoflurane. The customer may call 1-800-288-2121 for details.

If you have any questions regarding this COMHUB upgrade, contact Jim Wetzel at extension #2325.